

Pediatric Tracheostomy and Use of the Passy Muir® Valve

Course Outline

- Pediatric Tracheostomy
- Clinical Complications of Tracheostomy
- Clinical Benefits
- Assessment Guidelines
- Treatment Principles
- Therapy Techniques- with video examples

Pediatric Tracheostomy

- Timing of Pediatric Tracheostomy
 - Majority are trached under 1 year of age
 - Average age 2-3 years old
 - Critical time for development
 - Practice and hear their own vocalizations
- Indications for Pediatric Tracheostomy
 - Long term ventilation
 - Neuromuscular
 - Ventilator dependency
 - Upper airway obstruction
 - Stenosis
 - Laryngomalacia, tracheomalacia, bronchomalacia
 - Craniofacial malformations and syndromes
 - Vocal cord paralysis
- Pediatric Airway Anatomy
 - Softer, larger palate
 - Lumen size is smaller
 - Larynx is higher
 - Close approximation of larynx and tongue
 - Less cartilage (softer)
- Normal Trachea
- Grade 3 Subglottic Stenosis
- Grade 4 Complete Obstruction
- Laryngomalacia
- Vocal Cord Paralysis

Clinical Complications of Tracheostomy

- Effects of Tracheostomy on Communication Development
 - Caregiver interaction
 - Voice quality
 - Speech
 - 61% have articulation error

- Consonant and vowel errors
 - Excessive use of phonological processes (stridency and liquid deviation, cluster reduction, fronting)
- Language
- Impact on Swallow
 - Laryngeal tethering
 - Decreased sensation in the oropharynx
 - Reduced airway closure
 - Reduced subglottic air pressure
- Pediatric Studies
 - Swallow Physiology
 - 91% of 36 trach infants had swallowing disorders; half attributed to underlying neurological or anatomic deficits
 - Higher incidence of enteral feeding
 - Trach effect (Abraham (2000))
 - Toddlers with trach
 - Reduced laryngeal elevation
 - Delayed swallow response penetration
- Decreased Secretion Control and PEEP
 - Removal of filtration & humidification
 - Decreased cough
 - Abraham (2009)
 - All 50 children in study had secretion management issues
 - 98% level of trachea
 - 40% level of larynx
 - 56% oral
- Restore Airflow to Upper Airway
- Clinical Benefits
- Improved Secretion Management
 - Abraham (2009)
 - 24/49 children wearing Passy Muir valve full time
 - Secretion management within normal limits in average of two weeks
- Cotreatment
- Expedites Decannulation and Weaning
 - Restoration of normal physiology
 - Utilization of expiratory muscles
 - Accustomed to more normal breathing pattern
 - Increased PEEP and oxygenation
 - Able to communicate
 - Less WOB, easier to tolerate
- Improved Quality of Life

Types of Passy Muir Valves

- PMV® 2000 (clear) & PMV 2001 (Purple Color™)
- PMV Secure It®
- PMA 2000 Oxygen Adapter
- PMV 007 (Aqua Color™)
- Patient Care Kit
- Team Approach
 - Multi-disciplinary
 - Patient centered
 - Policies and procedures
- General Differences
 - Psychological aspects
 - Habilitation
 - Behavioral considerations
- Patient Selection
 - Cognitive status- awake, responsive, attempting to communicate
 - Medically stable
 - Able to tolerate cuff deflation
 - Able to manage secretions
- Factors Affecting Upper Airway Patency
 - Size of tracheostomy tube
 - Presence and degree of obstruction
 - Edema
 - Secretions
- To Assess for Upper Airway Patency
 - Deflate cuff
 - Finger occlude and voice or cough on exhalation
 - Use mirrors, cotton, feathers, whistles or bubbles to assist with the oral exhalation process
- Placement Guidelines
 - Patient education
 - Patient position
 - Suctioning
 - Achieve complete cuff deflation
 - Use the warning label provided with packaging
- Baseline Measurements
 - Oxygenation
 - Vital Signs
 - Breath Sounds
 - Color
 - WOB
 - Patient Responsiveness
- Distress... Signs and Symptoms

- Increased RR, HR
- Anxiety and fear
- Restlessness
- Increased irritability
- Stridor
- Grunting (infants)
- Retractions
- Nasal flaring
- Head bobbing
- Sniffing position
- Decreased BS during auscultation
- Decreased chest movement
- Decreased LOC
- Decreased PaO₂ (SaO₂)
- Increased PaCO₂
- Paleness or cyanosis
- Decreased perfusion/mottling
- Bradycardia/hypotension (this is a late sign)
- Guidelines
 - Approach to Education
 - Higher incidence of airway obstruction
 - Normal baselines
 - Reaction time
 - Transitions may be gradual
- Children 0-6 months old
 - Evaluation:
 - Infant communication developmental scale
 - Oral-motor function and feeding skills
- Treatment:
 - Parent education
 - Early developmental milestones
 - Non-nutritive oral stimulation/sucking
 - Bottle feeding
- Children 6-24 months old
 - Evaluation:
 - Receptive and expressive language
 - Play and cognition
 - Oral feeding Skills
 - Treatment:
 - Parent education
 - Facilitate oral exhalation
 - Facilitate vocalization
 - Augment with alternative communication systems as necessary

- Oral feeding
- Play, Play, Play!
- Goal #1 Increase Oral Exhalation
 - Techniques
 - Imitation
 - Blowing
 - Bubbles
 - Whistles
 - Horns, kazoos
 - Pinwheels
 - Straws
 - Cotton balls
- Toby Tracheasaurus™
- Goal #2 Increase Voicing
 - Activities
 - Planes, trains, and automobiles
 - Play dough
 - Rice and beans
 - Animals
 - Books
 - Songs
- Children 3 years old-School Age
 - Evaluation:
 - Receptive and expressive language
 - Oral feeding Skills
 - Treatment:
 - Parent education
 - Coordinate services with school, IEP
 - Elimination of negative behaviors
 - Intensive voice and speech therapy
- Additional Voice and Speech Goals
 - Improve voice quality
 - Weak
 - Harsh, hyper-functional
 - Pitch
 - Improve coordination of respiration and speech
 - Improve articulation and phonology
- Swallowing and Feeding Goals
 - Passy Muir valve for assessment (clinical & instrumental)
 - Passy Muir valve for oral feeding
 - Decrease sensitivity (facial/oral stimulation)
 - Tastes

- Textures
 - Feeding Equipment
 - Positioning
- Care, Cleaning and Lifetime of the Passy Muir Speaking Valves
 - Average lifetime of two months

Educational Opportunities

WEBINARS or SELF STUDY COURSES

Application of Passy Muir Swallowing and Speaking Valves Interdisciplinary

Tracheostomy Team: Where Do I Start?

Ventilator Basics for the Non-RT

Ventilator Application of the Passy Muir Valve

Pediatric Tracheostomy and Use of the Passy Muir Valve Pediatric Ventilator

Application of Passy Muir Valve

Swallow Function: Passy Muir Valve Use for Evaluation & Rehabilitation

www.PassyMuir.com

Passy-Muir, Inc. is an approved provider of continuing education through ASHA, AARC, and the California Board of Nursing.