

## Progression to Tracheostomy Decannulation: Role of the Speech-Language Pathologist

Suzanne E. Johnston, MA, CCC/SLP  
Kristin A. King, PhD, CCC/SLP



## Financial / Non-Financial Disclosures



- ◆ **Suzanne Johnston, MA, CCC-SLP**
- ◆ Financial: Part-time, clinical consultant for Passy Muir, Inc.
- ◆ No relevant non-financial disclosures



- ◆ **Kristin A. King, PhD, CCC-SLP**
- ◆ Financial: Full-time, Vice President of Clinical Education and Research for Passy Muir, Inc.
- ◆ No relevant non-financial disclosures

Disclosure: This presentation will focus primarily on the no-leak Passy-Muir® Valve and will include little to no information on other speaking valves.

## Learning Objectives

- ◆ Identify indicators for tracheostomy
- ◆ Understand anatomical and physiological consequences of tracheostomy
- ◆ Understand unique and intersecting roles of tracheostomy team members
- ◆ Understand the Scope of Practice and role of the SLP in identification and management of patient needs regarding decannulation
- ◆ Identify criteria for decannulation

## Indications for Tracheostomy

- ◆ Prolonged mechanical ventilation
- ◆ Inability to perform trans-laryngeal intubation
- ◆ Upper airway obstruction
- ◆ Secretion management
- ◆ Neuromuscular disease
- ◆ Respiratory compromise: ARDS, COPD



## Reputed Benefits of Tracheostomy

- ◆ Improved patient comfort/less need for sedation
- ◆ Lower WOB/faster weaning from MV
- ◆ Improved safety
- ◆ Improved oral hygiene and oral intake
- ◆ Less long term laryngeal damage
- ◆ Lower VAP rates
- ◆ Lower mortality
- ◆ Reduced ICU and overall LOS
- ◆ Earlier ability to speak/improved participation


Durbin, C., (2010). Tracheostomy: Why, when, and how? *Respiratory Care*, 55(8):1056.

## Metal Tracheostomy




## Air Filled Cuffs

- ◆ Cuff inflated
- ◆ Cuff deflated




## Water Filled Cuffs

- ◆ Cuff inflated
- ◆ Cuff deflated



## FOME-Cuff®

Self sealing



CONTRAINDICATED For Passy-Muir Valve USE

## Tracheostomy tubes (sizing)

BRAND	#	ID (MM)	OD (MM)	LENGTH (MM)
SHILEY	6	6.4	10.8	76
	8	7.6	12.2	81
PORTEX	6	6	8.3	55
	8	8	11	76
BIVONA	6	6	8.7	100
	8	8	11	120


- ◆ Consensus: the trach tube should be no larger than 2/3 of the inner diameter of the tracheal lumen (*adult standard*)

## Policies and Protocols

- ◆ Decannulation/weaning efforts begin *at time of intubation*
- ◆ Establish:
  - ◆ Effective humidification
  - ◆ Patency of airway
  - ◆ Appropriate cuff pressures and cuff deflation
  - ◆ Effective secretion management
  - ◆ Effective ventilation
  - ◆ Weaning plan and coordination of care in place for short term, as well as long term, needs


## Impact on Anatomy & Physiology

- ◆ Airway/ventilation
  - ◆ Separation of upper & lower airway
  - ◆ Change in pressures
    - ◆ PEEP, affecting gas exchange
    - ◆ Subglottic pressure
    - ◆ Oral, pharyngeal, esophageal pressures
  - ◆ Disuse muscle wasting/atrophy
  - ◆ Secretion management
  - ◆ Sensory changes
- ◆ Cuff complications



## Impact on Anatomy & Physiology

- Communication
  - Little to no voicing
  - Poor breath support
- Swallowing changes
  - Reduced subglottic pressure
  - Potential laryngeal tethering
  - Sensory awareness



## Trach Problems Identified by SLPs

- Lack of standardization for:
  - Sizing of the trach tube
  - Downsizing
  - Decannulation
  - Cuff deflation
- Trach care
- Use of Passy Muir Valves
- Referrals for swallowing assessments
- Limited patient/family education and discharge teaching




This Photo by Unknown Author is licensed under CC BY-SA

## Challenges to Establishing Teams

- No clinical consensus for decannulation protocols
- Members may vary per facility
- Defining team members' roles based on scope of practice and facility requirements
- Communication between members
- Establishing criteria for decannulation
- Coordination/timing of treatment of care


## WHO? Team Management of Tracheostomized Patients

- Not All Teams are Created Equal!



## Research: Supporting the Team

- Tobin, AE and Santamaria JD. (2008) An Intensivist-led Trach Review Team is Associated with Shorter Decannulation Time and Length of Stay: a prospective cohort study. *Critical Care*. 12 (2):R48. doi: 10.1186/cc6864
- Zaneta et al. (2014). Tracheal decannulation protocol in patients affected by Traumatic Brain Injury. *International Archives of Otorhinolaryngology*, 18(2): 108–114.
  - Faster decannulation, fewer weaning attempts
  - Reduced LOS
  - Cost savings
- Garrubba et al. (2009). Multi-disciplinary care for the tracheostomy patients: a systematic review. *Critical Care*. 13:R177
  - Reduction in time to decannulation
  - Shorter LOS
  - Fewer adverse events



## Defining SLP Role

- “The role of the SLTs [SLPs] is key, not only in assessing and managing swallowing and communication needs, but in contributing experience and expertise to all relevant tracheostomy-related decisions, as part of the MDT process.”
- McGrath (2014) The UK National Tracheostomy Safety Project

## ASHA Scope of Practice

- ◆ “Each practitioner evaluates his or her own experiences with preservice education, practice, mentorship and supervision, and continuing professional development. As a whole, these experiences define the scope of competence for each individual. The SLP should engage in only those aspects of the profession that are within her or his professional competence.”
- ◆ ASHA Scope of Practice for Speech Language Pathology, 2016
- ◆ Responsible for:
  - ◆ Optimizing a patient's ability to communicate and swallow, thereby improving QOL
  - ◆ Decisions are based on best available evidence
  - ◆ Work collaboratively

## WHAT? SLP Role in Decannulation

- ◆ Role in relationship to other team members
  - ◆ Create role-specific protocols for EACH member of the team
  - ◆ Consider establishing regular treatment times for streamlined scheduling
- ◆ Consistent education of all team members, on all shifts, who will be providing care

## SLP Role with Patients with Tracheostomy

- ◆ Identify established decannulation indicators pertinent to SLP Scope of Practice:
  - ◆ Level of alertness/cognitive and emotional state
  - ◆ Assessing for patent upper airway
  - ◆ Tolerance of cuff deflation
  - ◆ Use of speaking valve
  - ◆ Evaluating swallowing and secretion management
  - ◆ Identify and define aspiration risk



## Cuff Deflation Benefits

- ◆ Reduces aspiration<sup>1,2</sup>
- ◆ Improves laryngeal elevation
- ◆ Weaning time shorter with cuff deflation -avg of 3 days vs 8 days<sup>3</sup>
- ◆ Fewer respiratory infections, including VAP in cuff deflated group (20% vs. 36%)<sup>3</sup>
- ◆ Swallowing better in cuff deflated group and improved more from baseline<sup>3</sup>



1. Davis, et al. (2002). Journal of Intensive Care Medicine. 17(3): 132-135.  
 2. Ding, R. & Logeman, J. (2005). Head & Neck. 27(9):809-13  
 3. Hernandez, et al. (2013). Intensive Care Medicine. 39(6):1063-70

## How Does The Valve Work?

- ◆ Patented “no leak” design
- ◆ Opens only during active inspiration
- ◆ Closes at end inspiration
- ◆ Remains closed t/o expiratory cycle
- ◆ Air is re-directed thru the upper airway
- ◆ Offers a buffer to secretions



## Passy Muir® Valve Improves Weaning and Decannulation

- ◆ Improved scores on PAS<sup>1</sup>
- ◆ Restores expiratory airflow<sup>2</sup>
- ◆ Improves laryngeal clearance<sup>2</sup>
- ◆ Improved secretion rating scale<sup>3</sup>
- ◆ Maintains lung volumes<sup>4</sup>
- ◆ Restores subglottic pressure for cough<sup>5</sup>
- ◆ **Decreased Decannulation time**



1. Suiter, D. Head and Neck. 2005. Sep;27(9):809-13  
 2. Piguenti, Helene. Intensive Care Med. 2012 Jun;38(1):85-90.  
 3. Blumentfeld, L. Oral Abstract Presented at DRS Annual Meeting 2012  
 4. Gross, R., et al. (2006). The Laryngoscope, 116:753-761  
 5. Eibling, D. & Gross, R. (1996). Annals of Otolaryngology, & Laryngology, 105(4):253-8.




## WHEN? Predictors of Success

- ◆ Literature supports that there ARE established predictors of decannulation success:
  - ◆ Alert and responsive/level of consciousness
  - ◆ Resolution of tracheostomy indication
    - ◆ No acute respiratory compromise
    - ◆ Vent settings support weaning
  - ◆ Medically/hemodynamically stable

## WHEN? Predictors of Success

- ◆ Patent upper airway
- ◆ Tolerates cuff deflation, speaking valve, capping or plugging
- ◆ Good secretion management
  - ◆ Effective expectoration
  - ◆ Protective reflexes: cough and throat clear
  - ◆ Improved swallow to prevent aspiration
- ◆ Supportive environment post-decannulation

## SLPs and Treatment: RMST



## Sebastian



## Evaluation for Decannulation


- ◆ Reason for tracheotomy has resolved
- ◆ Medically stable
- ◆ Patent upper airway
- ◆ Tolerates speaking valve
- ◆ Can manage oral and tracheal secretions
- ◆ Tolerates capping/plugging
- ◆ Risk of aspiration assessed

## Candidates for Decannulation

- ◆ Consider original reason for the trach
- ◆ Weaned from mechanical ventilation, effective cough, no significant upper airway lesion<sup>1</sup>
- ◆ Absence of distress, stable arterial blood gases, hemodynamic stability, absent fever<sup>1</sup>
- ◆ A peak cough flow of 160 liters/minute<sup>2</sup>
- ◆ Survey: patient's level of consciousness, cough effectiveness, secretions, oxygenation<sup>3</sup>



1. Christopher, K. (2005). Respiratory Therapy. 50(4):538 -54.  
 2. Bach & Saporito, (1998). Chest. 110(6): 1566-71.  
 3. Sciflow, H. et al (2009). Respiratory Care. 54(12): 1658-68.



## Outcome Measures by SLP

- ◆ Establish patient's ability to:
  - ◆ Tolerate cuff deflation
  - ◆ Use upper airway for respiration, cough, throat clear, and speech
  - ◆ Manage secretions
  - ◆ Exhibit voice: quality
  - ◆ Swallowing safely
  - ◆ Participate in care
  - ◆ Understand education



## Factors for Decannulation Success

- ◆ Factors affecting weaning
  - ◆ Patient status
  - ◆ Tube size
  - ◆ Need and use of cuff
- ◆ For decannulation success
  - ◆ Sufficient air movement through upper airway when cuff is deflated or with uncuffed tracheostomy tube
  - ◆ Cuff must be *completely* deflated; open fenestration is insufficient for adequate air movement



## Air Whoosh - Backpressure



## Input from Team Members for Decannulation


- ◆ RT- secretion status, cough ability, airway patency, respiratory condition
- ◆ SLP- secretion management, cough, airway patency, swallow status
- ◆ RN- secretion status, level of consciousness
- ◆ MD- medical stability, whole patient assessment, recommendations from team

**All members communicate to each other regarding the patient's readiness to decannulate**

## Decannulation: Removal of Tracheostomy Tube



## Why a Trach Team: Summary



1. Communication
2. Patient Safety
3. Risk of Aspiration
4. Risk Associated with Trach Tube
5. Infection Control
6. Mechanical Ventilation
7. Long-Term Trach Placement
8. Education
9. Staff Confidence/Knowledge
10. Plan of Care and Continuity of Care
11. Quality of Care
12. Quality of Life

## Questions/Comments?

- ◆ What are your experiences – successes and frustrations, solutions with Trach Team work?
- ◆ Advice that you would share regarding your work with tracheostomized/ventilated patients?
- ◆ Recent research/developments you might share regarding SLP role regarding decannulation?

## RESOURCES



BIBLIOGRAPHY

<http://passymuir.com/asha2017>

## References

- ◆ Available at:  
[www.passymuir.com/ASHA2017](http://www.passymuir.com/ASHA2017)