Modes of Ventilation – Method or way a breath is delivered.

- **Control mode** Every breath delivered to patient is a mechanical breath. The breath may be triggered by a timing mechanism or patient effort. May be volume or pressure controlled. Examples: AC, VC/IMV, PC/IMV, PRVC.
- Supported or spontaneous mode Every breath is spontaneous, patient triggered, and supported by ventilator. Example: PS, CPAP
- **Combination mode**. Combination of both controlled and supported breaths. Example SIMV/PS is a combination of a set number of volume controlled breaths and pressure supported breaths.
- VC/AC Volume Control/Assist Control the ventilator is set to deliver a specific volume each breath, regardless of the amount of pressure required to deliver the volume. The clinician can set a high-pressure limit.
- SIMV Synchronized Intermittent Mandatory Ventilation tidal volume and rate are set, but the ventilator senses patient effort and "reschedules" mandatory (set) breaths.
- PS **Pressure Support** is a patient-initiated breathing mode in which the ventilator supports patient effort. Provides a set amount of pressure during inspiration to help patient draw in a spontaneous breath. The patient regulates the breath rate. The inspiratory time and/or volume of each breath may vary.
- PC **Pressure control** is an alternative to volume control. A pressure level is preset. Breaths are delivered at a preset frequency rate and inspiratory time. Pressure is constant throughout the delivered breath. Tidal volume can be variable. Patient can breathe above set rate.
- PEEP Positive End Expiratory Pressure maintains small end-expiratory pressure to help prevent alveolar collapse and improve oxygenation. Patients are often started on 5 cmH2O of PEEP.
 CPAP Continuous Positive Airway Pressure is positive pressure maintained in the airway to prevent alveolar collapse, but is generally used to describe positive pressure set in spontaneous mode.

<u>Settings</u> – in addition to the mode of ventilation, the following are physician ordered and/or set by the clinician.

- V_T- **Tidal Volume** Volume of air delivered per breath. Pressure to deliver the breath may vary depending on lung mechanics.
- RR/F **Respiratory Rate/Frequency** Set frequency of ventilator delivered breaths per minute.
- I-Time **Inspiratory Time**, expressed in seconds, is the amount of time spent in inspiration during the total respiratory cycle. Example: RR is 12, total cycle time is 5 seconds (60/12 = 5 seconds). If I-time set at 1 second, then 1 second is spent in inspiration, 4 seconds are available for exhalation.
- FIO2 Fraction of Inspired Oxygen. the amount of oxygen the ventilator delivers, expressed as a percentage. Room air is 21%.
- Sensitivity The level of effort from the patient needed to "trigger" the ventilator to deliver a breath from the ventilator. Increase sensitivity to decrease patient effort.

<u>Measured</u> – The following parameters are measured by the ventilator and displayed on the ventilator monitor.

Exhaled Tidal Volume - amount of air exhaled each breath

Exhaled Minute Volume - amount of air exhaled in one minute

Peak Inspiratory Pressure – maximum amount of pressure reached during inspiration

RR/F - Total respiratory rate includes set rate and patient initiated breaths

<u>Alarm Settings</u> –

High and low pressure High and low exhaled tidal volumes High respiratory rate