



GdYU_]b['J U`j Y`I gY`k]hΛ`
 HfUWYcghca mUbX`
 A YWUUb]WU`'J Ybh]`Uh]cb.`
 Bck ž`K YUY`HU`]b[°


DfYgYbhYX`]b`&\$&)

© 202) Passy-Muir, Inc.

These materials are protected by the Copyright laws of the United States and may not be copied, distributed, or republished by the recipient. Recipients are licensed to use the materials for their personal use as an attendee of this seminar.

Passy-Muir®, PMV®, the purple valve®, and the aqua valve™ are trademarks of Passy-Muir, Inc.

Speaking Valve Use with
Tracheostomy and Mechanical Ventilation:
Now, We're Talking!




1

Welcome!

- In-Person Considerations
 - Social awareness
 - Networking options
 - Presentations
 - Groups: Purple and Aqua
 - Hands-On
 - Demonstrations
 - Case studies/mock assessments
 - Panel discussions
 - Food and drink

- Chat with the clinicians!
- Meet fellow professionals!
- Learn about new resources!



2

Speaker Disclosures:

Financial: Full-time with Passy-Muir, Inc.
Non-financial: No relevant non-financial disclosures

Kristin A. King, PhD, CCC-SLP
Vice President of Clinical Education and Research



Gabriela Ortiz, BSRT, RCP
Clinical Specialist



3


Disclosure

Portions of the information in this seminar relate to products of Passy-Muir, Inc. The content of the seminar is for your general educational information only. Information you learn cannot replace your clinical judgment or the judgment of treating physicians and other healthcare professionals working with patients. The information provided in this seminar should not be considered as medical advice.

4

Seminar Overview

Fundamentals: From Tracheostomy Tubes to Mechanical Ventilation Terminology



- 7.5 hours today

- Up to 6 recorded hours
 - 1 pre-requisite webinar hour
 - 5 post-seminar recorded webinar hours

- 13.5 total hours, if all are completed.


5

Pre-requisite webinar

- *Fundamentals: From Tracheostomy Tubes to Mechanical Ventilation Terminology*
- Overview of:
 - Tracheostomy tubes
 - Cuff types and management
 - Considerations as it relates to speaking valve use
 - Review of terminology related to mechanical ventilation
 - Implications for patient management

6

Seminar Participation



- Interactive questions
- Q & A
- Now, get your phones out.
- Let's try it!

7

slido




What is the actual color of the PMV 007?


Start presenting to display the poll results on this slide.

8

Assessment and Placement: Non-Ventilator




9



10

Begins in the ICU: Effects of Bed Rest




- The negative impact of bed rest is well known
- No evidence supports efficacy of bed rest
- Disuse atrophy at the cellular level begins within 4 hours of implementing bed rest
- Healthy adults, bed rest¹
 - Strength declined by 1 – 1.5% per day
 - Mood changes
 - Loss of coordination, balance and work tolerance
 - Casting: Strength declines by 25% in 7 day²

• Griffiths et al. Nutrition 1996; 11:428-432
• De Jonghe et al. ICCM 2000; S309-315

11

Endotracheal Tube Considerations


- Intubation
 - Why?
 - Respiratory compromise due to injury or illness
 - Potential for respiratory compromise
 - Potential side effects
 - Mucosal injury
 - Injury to the vocal folds
 - Tracheal edema, ulceration, and stenosis



12

Physiologic Changes after Tracheostomy


Respiration: Patient inhales and exhales through open trach tube. No airflow past inflated cuff



13

Physiologic Changes after Tracheostomy


- Speech
- Smell
- Taste
- Sensation
- Reduced positive airway pressure
 - Poor secretion management
 - Reduced cough



14


Cuff: Choices and Management

- Cuff up or down?
 - Purpose of cuff
 - Cuffs and aspiration
- Cuff pressures
 - 20 - 25 cmH₂O
 - Minimal leak
 - Minimal occlusion



15

Reminder: Passy-Muir Valve



Bias-closed


No-leak

PMV 2001 (Purple color™)

What are the benefits of the Valve?

19

slido




What is a benefit of using the Passy Muir Valve?

① Start presenting to display the poll results on this slide.

20

Benefits of Closing the System with PMV

- Restores normal physiology
- Reconnects the upper and lower airway
- Providing a closed system
 - Communication
 - Smell and taste
 - Secretion management
 - Sensation
 - Cough
 - Swallowing
 - Positive airway pressure
 - Quality of life



21

Psychological Benefits

Communication with family

Participation in decision making

Reduced sense of isolation/ anxiety

Better sense of well-being


Communication with caregivers

• Freeman-Sanderson, A. L., Tagher, L., Elms, M. R., & Kenny, B. (2018). Quality of life improves for tracheostomy patients with return of voice: A mixed methods evaluation of the patient experience across the care continuum. *Intensive Critical Care Nursing*, 46, 10-16. doi:10.1016/j.iccn.2018.02.004
• Freeman-Sanderson, A. L., Tagher, L., Elms, M. R., & Phipps, P. R. (2016). An intervention to allow early speech in ventilated tracheostomy patients in an Australian intensive care unit (2019): A randomized controlled trial. *Australian Critical Care*, 30(3), 134. doi:10.1016/j.aucc.2015.12.012
• Freeman-Sanderson, A. L., Tagher, L., Elms, M. R., & Phipps, P. R. (2016). Quality of life improves with return of voice in tracheostomy patients in intensive care: An observational study. *Journal of Critical Care*, 31, 198-204. doi:10.1016/j.jcc.2016.01.012
• Freeman-Sanderson, A. L., Tagher, L., Elms, M. R., & Phipps, P. R. (2016). Return of voice for ventilated tracheostomy patients in ICU: A randomized, controlled trial of early-targeted intervention. *Critical Care Medicine*, 44(5), 1079-1081. doi:10.1097/ccm.0000000000001010

22

Impact on PEEP


- Closed System vs Open
 - Improved gas exchange
 - Improved oxygen saturation levels
 - Decreased risk of atelectasis
- “My patient is not tolerating cuff deflation trials”



23

Initiating the Assessment: Team Approach


- Have a plan: Who does what
- Block off the time
- Education
- Reassure the patient
- Perform good oral care
- Suctioning as needed
- Body position and posture
- Position of head, neck, and tracheostomy tube



24

Patient Selection


- Awake and alert
- Medically stable
- Complete cuff deflation
- Manageable secretions
- Patent upper airway



25


Checklist: Take Baseline Measurements

- Oxygenation
- Vital signs
- Breath sounds
- Color
- Work of breathing
- Patient responsiveness




26


Education




27



28




29



30


Deflate Cuff



31

Assess Upper Airway Patency


- Deflate cuff
- Ask patient to inhale
- Finger occlude and speak or cough on exhalation
- Transtracheal pressure measurements




32

Assessment for Placement


- Transtracheal pressure measurements
 - Back pressure
 - Air trapping
 - Assessing for patent upper airway




33



34




35



36


Initial Placement: Voicing



37

Advantages of a Closed Respiratory System vs Open Tracheostomy


- Open tracheostomy
 - Reduced airflow
 - Reduced positive airway pressure
 - Reduction in the pressurized system
- Closed Respiratory System
 - Allows graded exhalation and pressure regulation
 - Feeding and Swallowing
 - Posture and balance
 - Upper extremity force/strength



38

Pressurized System


- Restored or improved pressurized system:
 - Intraoral
 - Subglottic pressure
 - Respiratory – PEEP
 - Esophageal ??
 - Intrathoracic
 - Respiratory
- Leads to improved:
 - Feeding and swallowing
 - Cough and throat clear
 - Trunk support and postural control
 - Respiratory function



39


Wear Time

- Patient specific
 - Patient's cognitive status
 - Medical needs
- Minutes to hours
- Treatment plan



40


More Education





41

Care and Cleaning

- Average lifetime of 2 months









42

Decannulation



43


Pre and Post Decannulation



44

In Summary: What has been accomplished?

- Early intervention:
 - Avoid disuse atrophy
- Close the system to improve:
 - Phonation: access to vocal communication
 - Sensation and secretion management: cough and throat clear
 - Taste and smell
 - Time to weaning and decannulation



45

Troubleshooting and Treatment

46

Cuff issue?

47

Transitioning and Troubleshooting

What if the patient cannot voice or has difficulty?

What are the steps to consider?

Step 1: Check the cuff status and patient positioning


Step 2: Reassess airway patency

Step 3: Assess the stoma/leak

Step 4: What to do next?


48

No Voice

A photograph of a man lying in a hospital bed. He is wearing a patterned hospital gown and has a Passy-Muir tracheostomy tube inserted into his tracheostomy. He is looking towards the camera with a slight smile. Medical equipment and a blue chair are visible in the background.


49

Too Many Secretsions

A photograph showing a healthcare worker wearing a white protective gown, gloves, and a face mask. They are leaning over a patient who is lying in a hospital bed. The patient appears to be wearing a tracheostomy tube. The background shows a typical hospital room setting with medical equipment.


50

Case Study: Troubleshooting

A photograph of two healthcare workers, a man and a woman, both wearing blue gloves, attending to a patient lying in a hospital bed. The patient is wearing a white hospital gown and has a tracheostomy tube. The man is on the left, and the woman is on the right. They appear to be adjusting the patient's position or the tube. Medical equipment is visible in the background.

51

slido




What were some signs of distress that the patient was showing?

① Start presenting to display the poll results on this slide.

52

slido




What are possible causes of back pressure?

① Start presenting to display the poll results on this slide.

53

Factors Affecting Expiratory Air Flow

- Size or type of tracheostomy tube
- Presence and degree of obstruction
- Edema
- Secretions
- Incomplete cuff deflation
- Tube position



54

Airway Obstruction



55

Troubleshooting:
Downsize or Different Brand Tube

Trach A

Trach B



56

Case Study: Resolution



57

Questions to Determine Therapy

- What is diagnosis?
- Why do they have difficulty with:
 - Voice?
 - Breath support?
 - Language and/or cognition?
 - Dysarthria?
- What about swallowing?


58

Goals

- Wear time
 - Patient will wear the speaking Valve:
 - For ____ minutes to improve communication.
 - During ____-minute therapy session without need for Valve removal.
 - For ____-hour periods of time while awake.
- Other goal areas that impact wear time
 - Participate in conversation with audible voicing on ____ out of ____ sentences.
 - Complete ____ number of RMT tasks while wearing the speaking Valve.


59

Break: 15 Minutes



60


Basics of Ventilator Application
of the No-leak Valve



61

Indications for Invasive Mechanical Ventilation

- Can no longer support with NIV
- Airway protection
- Hypercapnic respiratory failure
- Hypoxemic respiratory failure
- Cardiovascular distress
- Anticipated patient decline or impending transfer





62

Positive Pressure Ventilation

- The primary approach to mechanical ventilation is positive pressure
- The flow's driving pressure is higher than the pressure inside the lungs
- The pressure is high enough to overcome resistance

Negative Pressure vs Positive Pressure Ventilation






63

Invasive Ventilation


- Usually requires airway to be sealed with little to no leak present.
- Seal is achieved with a cuff at the end of the artificial airway.




64

Invasive Ventilation

- Endotracheal Tube




- Tracheostomy Tube



65

Modes of Ventilation

- Control Modes
 - Continuous Mandatory Ventilation (CMV)
 - Assist Control (A/C)
- Spontaneous Modes
 - Pressure Support (PS)
 - Continuous Positive Airway Pressure (CPAP)
- Combination Mode
 - Synchronized Intermittent Mandatory Ventilation (SIMV)



66

Ventilator Settings


Requiring physician order:

- Mode (A/C, SIMV)
- Breath Type:
 - Pressure Control (PC)
 - Volume Control(VC)
 - Pressure Support (PS)
 - Pressure Regulate Volume Control(PRVC)
- Frequency(f)/Respiratory Rate(RR)
- Positive End-Expiratory Pressure (PEEP)

Therapist driven:

- FiO₂
- Flow rate
- Trigger sensitivity
- Inspiratory time/flow limits
- Alarm settings

67



Measured Ventilator Parameters

- Exhaled Tidal Volume
- Exhaled Minute Volume
- Peak Inspiratory Pressure (PIP)
- Total Respiratory Rate
- Mean Airway Pressure
- PEEP


68

Steps For In-line Valve Placement


69

Why use a Passy-Muir Valve with patients who are mechanically ventilated?

- Verbal communication
- Improved lung recruitment and diaphragm involvement
- More rapid weaning from the ventilator
 - Rehabilitation tool
- Improved secretion management
 - More effective cough
 - Reduces need for suctioning
- Improves quality of life



70




71

Step 1: Assessment

72

Patient Selection Criteria




- Awake and alert
- Hemodynamically stable
- Able to manage complete cuff deflation
- Manageable secretions
- Patent Airway

73

Establish Baseline:
Assess Vital Signs and Work of Breathing

- Oxygenation
- Vital Signs
- Breath sounds
- Color
- Work of breathing
- Patient responsiveness




74

Assess Ventilator Parameters


Three parameters that give you the general state of your patient's respiratory status:

1. FiO_2
2. PEEP
3. PIP



75

Assess Ventilator Parameters




FiO₂

- Fraction of inspired Oxygen
- Room Air 21%
- Supplemental O₂ > 21%

76

Assess Ventilator Parameters




PEEP

- Positive End-Expiratory Pressure
- Resistance to exhaled volume that creates back flow and pressure that stents the alveoli open
- PEEP and FiO₂ work together to improve oxygenation

77

Assess Ventilator Parameters




PIP (Peak Inspiratory Pressure)

- The max amount of pressure to deliver volume
- Sum of the inspiratory pressure required to deliver volume + PEEP
- PIP indicates the compliance of the lungs

78

Suggested Ventilator Parameters

- $\text{FiO}_2 \leq .50$
- $\text{PEEP} \leq 10 \text{ cmH}_2\text{O}$
- $\text{PIP} \leq 40 \text{ cmH}_2\text{O}$
- Any conventional mode of ventilation is compatible with the Valve.



The image shows a ventilator control panel with multiple waveforms (Pressure, Volume, Flow) and numerical readouts for parameters like PEEP (27), PIP (12), and FiO2 (3.5). The interface is dark with green and yellow highlights.


79

Step 2: Patient Preparation and Education

80

Team Approach

- Timing and tube selection
- Introducing a speaking valve
- When to downsize
- Plan of care
- Decannulation
- Impacts continuity of care
- Impacts safety, length of stay, and cost




A group of six healthcare professionals, including nurses and doctors, standing together in a professional setting. They are wearing scrubs and lab coats.

81

© Passy-Muir, Inc.

Patient Preparation


- Body position and posture.
- Position of head, neck, and tracheostomy tube.



82

Pre-Placement, General Observations, and other Considerations

- Have a plan and block time
 - Pick a good time of the day
 - Reduce noise and interference
- Education
 - Reassure the patient
- Address pain issues
- Position the patient




83

Step 3: Assess For Airway Patency

84

Airway Patency Assessment with Mechanical Ventilation


- Requires complete cuff deflation
- Assess the leak or airflow into the upper airway
- Use vent parameters to determine airway patency
- The type of breath matters



85

Cuff Deflation and Mechanical Ventilation


1. Set parameters do not change when cuff is deflated.
2. Cuff deflation generates less resistance to flow.
3. Ventilatory system is no longer sealed, there is a leak.



86

VC: Patient Assessment


PIP	25 cmH ₂ O
PEEP	5 cmH ₂ O
FiO ₂	30%
Vt	500 mL
Vte	500 mL



87

Ventilator Assessment


- Note Vent Settings:
 - Set Vt
 - PEEP
 - RR/f
 - FIO₂
- Note Vent Measurements
 - PIP
 - Exhaled Vt (Vte)
 - Total RR
 - MV



88

Upper Airway Patency Assessment

- Turn Down PEEP
 - PEEP down by 5
- Then,
 - Slow cuff Deflation




89

Upper Airway Patency

Cuff Inflated-Closed Circuit


PEEP	5 cm H ₂ O
FIO ₂	.30
PIP	25 cmH ₂ O
Vt	500 mL
Vte	500 mL



Cuff Deflated-Open Circuit


250 mL

PEEP	0 cmH ₂ O
FIO ₂	.30
PIP	10 cmH ₂ O
Vt	500 mL
Vte	250 mL




90

Upper Airway Patency Assessment



91


PC: Patient Assessment



92

Ventilator Assessment


- Note Vent Settings:
 - Set IP
 - PEEP
 - RR/f
 - FiO₂
- Note Vent Measurements
 - PIP
 - Exhaled Vt (Vte)
 - Total RR
 - MV



93

Upper Airway Patency Assessment

- Turn Down PEEP
 - PEEP down by 5
- Then,
 - Slow cuff Deflation




94

Upper Airway Patency


Cuff Inflated-Closed Circuit

Cuff Deflated-Open Circuit

PEEP	5 cmH ₂ O
FiO ₂	.30
PIP	25 cmH ₂ O
Vti	550 mL
Vte	550 mL



PEEP	0 cmH ₂ O
FiO ₂	.30
PIP	20 cmH ₂ O
Vti	700 mL
Vte	330




95


Step 4: Assemble the Necessary Parts & Pieces

96


© Passy-Muir, Inc.



97



98



99

Step 5: Place the Valve In-line and Assess the Patient

100

Vital Signs and Work of Breathing


- Oxygenation
- Vital Signs
- Breath sounds
- Color
- Work of breathing
- Patient responsiveness



101


Troubleshooting

- PIP increases with each breath
 - Increased WOB is observed
 - High pressure limit alarms
- Valve is removed
 - A loud whooshing sound heard from the tracheostomy tube
- What could cause this?
- Recommendations?



102


Listen for Back Pressure



103


Factors Affecting Airway Patency

- Tracheostomy tube
- Cuff issues
- Airway obstruction



104

Importance of Tube Position



105

FOME-Cuff®

Self-sealing

A yellow, bulbous, self-sealing device with a white connector and a red tube attached.


CONTRAINDICATED For Passy-Muir Valve USE

106

Step 6: Adjust the Vent as Necessary

107

Consider Switching to NIV


A white medical device with a screen displaying waveforms and a control panel with buttons and a dial.

- Approach to breath delivery is very similar to invasive ventilation.
- Difference is, there is almost always a leak present.

108

NIV


- Almost always a leak present.
- Flow loss makes it difficult to measure pressures and volumes.
- Most parameters are calculated instead of measured.



109

Considerations with NIV


- Airway patency assessment should not be done in NIV.
- Uses the same settings or as close as possible.
- Maybe necessary to do a trial prior to cuff deflation and Valve placement.



110

Ventilator Assessment and Adjustments

- Adjust PEEP
- Evaluate sensitivity
 - Pressure vs. Flow Trigger




111

Ventilator Assessment and Adjustments

• Volume compensation

• Increase V_T in small increments to achieve pre-cuff deflation PIP




112

Ventilator Assessment and Adjustments

• Pressure Ventilation

• May adjust to achieve audible voice and adequate ventilation



113

Ventilator Assessment and Adjustments

• Flow limit


• Increase the % flow deceleration

• Ranges 20 to 80%

• Time limit

• Set I-time

• 1 second for most adults



114

Alarm Settings – Safe Practice

Low exhaled Vt and Ve alarms

Low pressure alarm


- Set 5 to 10 cmH₂O below PIP

High pressure alarm

- Set 10 cmH₂O above PIP


High respiratory rate

- 10 to 15 above baseline




115

Ventilator Settings and Alarm Management



116

Case Study: Vent




117

© Passy-Muir, Inc.

Gil


- Ventilator settings:
 - A/C RR 8
 - V_T 700 mL
 - PIP 25 cmH₂O
 - PEEP 5 cmH₂O
 - F_iO₂ .28
- Tracheostomy
 - 1 month
 - Size 8 Shiley XLT



118

Cuff Deflation Assessment


- Adjust PEEP
- Slow cuff deflation
- Ventilator:
 - Exhaled V_T 300 mL
 - PIP 12 cmH₂O
- Patient:
 - Weak cough
 - Voicing
- Should the Valve be placed in-line?



119

Vent Changes Increase Success


- Valve is placed in-line
- Assessment reveals:
 - Whispers only
 - Poor chest expansion
 - Increased RR
- What ventilator change could be made?



120


Vent Changes Increase Success

- Vent change:
 - Increase V_T to meet but not exceed pre-cuff deflation PIP



121

Gil



122

LUNCH




& LEARN

123

Ventilator Demonstration

Ventilator Adjustments




PEEP	FiO ₂
5 cm H ₂ O	21
5 cm H ₂ O	21
5 cm H ₂ O	21
5 cm H ₂ O	21

PIP	V _T
20 cm H ₂ O	500 cc
20 cm H ₂ O	500 cc
20 cm H ₂ O	500 cc
20 cm H ₂ O	500 cc

124

Breakout Sessions:
Tracheostomy Tubes and PMVs, Cuff
Management, Vent Application, and
Mock Assessments



125

Placing a Valve In-Line



126

Tracheostomy Tubes and PMVs:
Purple Team




127

Cuff Management:
Aqua Team



128

Mock Assessments:
Purple and Aqua Teams




129

Case Study #1: Non-Vent

- 68 y/o male
- Admitted for complications related to OSA
- PMH: obstructive sleep apnea, hypertension, type II diabetes, hypercapnia, chronic pain, BMI > 40, peripheral neuropathy, GERD, nocturnal CPAP with full face mask, depression
- Surgical history:
 - Prior admission: Patient had tonsillectomy to improve airflow through the upper airway with full face CPAP mask.
 - Current: Surgical tracheotomy performed 4 days ago. Patient received #8 Shiley XLT, cuffed.
- Patient is awake and alert, communicating with finger occlusion and nonverbal communication (gestures and writing) during the day. Patient demonstrates frequent coughing and is suctioned regularly.
- Consult received for PMV assessment.

130

slido



What other information would you like to know?

① Start presenting to display the poll results on this slide.

131

slido

✓

-

○

-

Does this patient meet the criteria for PMV assessment?

① Start presenting to display the poll results on this slide.

132

slido




Put the following steps in order for this patient's assessment.

① Start presenting to display the poll results on this slide.

133

slido



What additional recommendations should be considered?

① Start presenting to display the poll results on this slide.


134

Case Study #2: Non-Vent


- 59 y/o male
- Admitted for surgical management of laryngotracheal stenosis
- PMH: laryngotracheal stenosis, lupus, HTN, GERD, depression, CAD, lupus nephritis
- Surgical history:
 - Laryngotracheal reconstruction
 - Suprastomal stent
 - Tracheotomy with #6 cuffless Shiley
- Patient is awake and alert, medically stable, demonstrating desire to communicate
- Consult received for PMV assessment

135

slido




What do you need to know more about to complete an assessment of this patient?


 Start presenting to display the poll results on this slide.

136

slido




What considerations would you have with a history of laryngotracheal reconstruction?


 Start presenting to display the poll results on this slide.

137

slido





What factors do you want more information about?

 Start presenting to display the poll results on this slide.

138

slido

Please download and install the Slido app on all computers you use



Does the position (suprastomal) of the stent make a difference?

⌚ Start presenting to display the poll results on this slide.


139

Case Study #3: Non-Vent

- 52 y/o female
- Admitted for SOB
- PMH: hypertension, type II diabetes, chronic kidney disease, SARS-CoV-2, dyspnea, airway stenosis
- Previous hospitalization: SARS-CoV-2, intubated for 26 days, discharged 1 week ago to inpatient rehabilitation
- Patient continued to complain of shortness of breath with rest and exertion, occasional stridor noted, crackling breath sounds
- Surgical history:
 - Laryngoscopy: airway stenosis
 - Surgical tracheotomy: # 6 Bivona TTS
- Chest X-ray: RLL pneumonia, aspirates
- MBSS: open tracheostomy tube, pharyngeal dysphagia with aspiration of thin liquids
- Supplemental oxygen at 3L via trach mask
- Consult received for PMV assessment.

140

slido




Do you have enough information to proceed with this assessment?

⌚ Start presenting to display the poll results on this slide.

141

slido




Put the following steps for this patient's assessment in the proper order.

① Start presenting to display the poll results on this slide.

142

slido



Which instrumental assessment would you consider for this patient?

① Start presenting to display the poll results on this slide.

143

Case Study #1: Vent

Patient history:


- Admitted to ICU s/p pneumonia w/ 2 failed extubations
- h/o smoking, HBP, mild COPD
- Tracheal 2 weeks earlier
- Shiley #7 cuffed

Clinical findings:

- Awake, alert, following commands
- HR and SpO₂ are normal
- Hemodynamically stable


Pre-Cuff Deflation:	Settings: AC/VC RR = 16 bpm Vt = 500 mL PEEP = 5 cmH ₂ O FIO ₂ = 40	Measurements: PIP = 25 cmH ₂ O RR = 18 bpm Vte = 500 mL
Post-Cuff Deflation:	Settings: AC/VC RR = 16 bpm Vt = 500 mL PEEP = 0 cmH ₂ O FIO ₂ = 40	Measurements: PIP = 17 cmH ₂ O RR = 18 bpm Vte =

144



Are the patient's vent parameters stable?

145



What should you adjust on the vent to avoid auto-triggering?

146

Case Study #2: Vent

Patient history:

- 62 y/o male. s/p intubation for COVID for 45 days; no significant PMH; bedside/perc trach 2 weeks ago.
- Admitted to step-down ICU 3 days earlier.
- Bivona #6, Fome-Cuff

Clinical plan and findings:

- Weaning from the vent
- Trial PMV
- SpO₂ = 96%
- Hemodynamically stable
- No s/s of respiratory distress

Pre-Cuff Deflation:


Settings:

AC/VC
RR = 16 bpm
Vt = 500 mL
PEEP = 5 cmH₂O
FIO₂ = .48

Measurements:


PIP = 25 cmH₂O
RR = 18 bpm
Vte = 499 mL

147



What vent adjustment is recommended before cuff deflation?

148




Will the PIP drop?

149

Trachlore, Barriers, and More

- Panel Discussion – starting point:
 - You have to wait until a patient is weaned from the ventilator.
 - Our patients are too sick to use a Valve.
 - You need a fenestrated tracheostomy tube.
 - We have to keep the cuff inflated due to aspiration.
 - My patient cannot tolerate cuff deflation trials, so they are not ready for a Valve.
 - My patient speaks with a leak, so a Valve is not needed.
- What have you heard?


150



Audience Q&A

151

Receiving CEUs for this Course



- You will have 5 days from the time this course ends to complete the evaluation, which is required to receive credit
- Go to: <https://ep.passy-muir.com>
- Login or create an account
- Click on the purple box
 - Upper righthand corner
 - Labeled "Enter Meeting Code Here"
- The meeting code is:

Enter Meeting Code Here

152
