



# **Clinical Insights for Assessing and Using No-Leak Speaking Valves with Trachs and Vents: Key Takeaways for Pediatric and Adult Patients**

**Learning Through Lecture, Case Studies, Mock  
Assessments, Hands-On Training, &  
Demonstrations**

DfYgYbhYX`]b`&\$&)

© 2021 Passy-Muir, Inc.

These materials are protected by the Copyright laws of the United States and may not be copied, distributed, or republished by the recipient. Recipients are licensed to use the materials for their personal use as an attendee of this seminar.

Passy-Muir®, PMV®, the purple valve®, and the aqua valve™ are trademarks of Passy-Muir, Inc.

## Clinical Insights for Assessing and Using No-Leak Speaking Valves with Trachs and Vents: Key Takeaways for Pediatric and Adult Patients



Learning Through Lecture, Case Studies, Mock Assessments,  
Hands-On Training, and Demonstrations

1

---

---

---

---

---

---

---

---

## Disclosures

Kristin King, PhD, CCC-SLP



Financial: Full-time Vice President  
of Clinical Education and  
Research, Passy-Muir, Inc.

Non-financial: Nothing relevant to  
disclose

Gabriela Ortiz, RCP, RRT



Financial: Full-time Clinical  
Specialist, Passy-Muir, Inc.

Non-financial: Nothing relevant  
to disclose

Passy-Muir, Inc., is a medical device manufacturer that produces and markets PassyMuir® Tracheostomy & Ventilator Swallowing and Speaking Valves for patients with ventilator dependence and non-ventilator tracheostomy. This presentation will focus primarily on the no-leak PassyMuir® Valve and will include little to no information on other speaking valves.

2

---

---

---

---

---

---

---

---

## Disclosure

The content of the seminar is for your general educational  
information only.

Information you learn cannot replace your clinical judgment or the  
judgment of treating physicians and other healthcare professionals  
working with patients. The information provided in this seminar  
should not be considered as medical advice.

3

---

---

---

---

---


---

---

---

## Seminar Overview

Fundamentals: From Tracheostomy Tubes to Mechanical Ventilation Terminology



- 7.5 hours today and 4 hours tomorrow = 11.5 hours total
- Up to 6 recorded hours
  - 1 pre-requisite webinar hour
  - 5 post-seminar recorded webinar hours
- 17.5 total hours, if all are completed.

4

---

---

---


---

---

---

---

## Seminar Participation



- Interactive questions
- Q & A
- Now, get your phones out.
- Let's try it!

5

---

---


---

---

---

---

---



## What is the actual color the PMV 007?

6

---

---

---

---

---

---

---

## Assessment and Placement of the Passy-Muir Valve: Non-Vent



7

---

---

---

---

---

---

---

## Indications for Tracheostomy

- Prolonged mechanical ventilation
  - Neuromuscular
  - Ventilator dependency
- Inability to perform trans-laryngeal intubation
  - Trauma, max/fax deformity
- Airway obstruction
  - Subglottic stenosis
  - Tracheal or laryngeal malacia
  - Vocal fold paralysis
- Secretion management

8

---

---

---

---

---

---

---

## Pre-trach



9

---

---

---

---

---

---

---

### When: Timing of Tracheostomy

- Number of days intubated
- Pre-planning
- Age of child



Does time of tracheotomy procedure have an impact on outcomes?

10

---

---

---

---

---

---

---

---

### Post-trach



11

---

---

---

---

---


---

---

---

### Vocal Tract Development

Vocal Tract Length (Birth)	Vocal Tract Length (Adulthood)
6-8 cm	14-18 cm
1:2 ratio	1:1 ratio



12

---

---

---

---

---


---

---

---

### Calculating Tube Size

- **ATS Consensus:**
  - Allows for airflow for speech
  - Colinear with trachea
  - 2 cm beyond stoma
  - 1-2 cm above carina
  - Peds are usually up-sized, not down-sized



13

---

---

---

---

---

---

---

### Calculating Tube Size

Age or Weight	Tracheostomy Tube Size
Babies who weigh less than 1000 g	2.5 mm
Babies who weigh 1000-2500 g	3 mm
Neonates aged 0-6 months	3-3.5 mm
Infants aged 6 months to 1 year	3.5-4 mm
Infants aged 1-2 years	4-4.5 mm
Children older than 2 years	[age in years + 16] divided by 4

14

---

---

---

---

---

---

---

### Neonate vs Pediatric vs Adult Trach Tube



15

---

---

---


---

---

---

---

### Flex-tend Trach Tube



16

---

---

---

---


---

---

---

---

### Airflow Changes with Tracheostomy



17

---

---

---

---

---

---

---

---

### Implications of Tracheostomy in NICU, Pediatrics, & Adults

- Loss of sense of taste and smell
- Reduced airway closure
- Altered subglottic pressure
- Decreased sensation
- Increased accumulation of secretions
- Altered awareness and management of secretions
- Altered cough

18

---

---

---

---

---

---

---

---

Millie



19

---

---

---

---

---

---

---

David



20

---

---

---

---

---

---

---

## Safety Issue

"Being unable to communicate is emotionally frightening for children and can lead to an increase in sentinel events, medical errors and extended lengths of stay"

• Costello, Patak, and Pritchard (2010). Communication vulnerable patients in the pediatric ICU: Enhancing care through augmentative and alternative communication. *Journal of Pediatric Rehabilitation Medicine: An Interdisciplinary Approach* 3: 289-301

21

---

---

---

---

---

---

---





22

---

---

---

---

---

---

---

---

### JCAHO Requires Communication

Meeting JCAHO's goal is not the only reason to encourage patients to have a voice in their care:

Those who are actively involved in their care often have better outcomes.

Joint Commission International Center for Patient Safety. (2006). Patient-inclusive care: Encouraging patients to be active participants in their care. Patient Safety Link, 2(2), 1. <http://ajvos.modernmedicine.com/modern-medicine/content/jcahos-patient-safety-goals-part-1-practical-guide?page=full>

23

---

---

---

---

---

---

---

---

### Application of the Passy Muir Valve



24

---

---

---

---


---

---

---

---

### Team Approach



25

---

---

---

---

---

---

---

### Open Tracheostomy Tube



26

---

---

---

---

---

---

---

### Passy-Muir Valves

- Bias-closed

- No-leak



PMV® 2001 (Purple color™)



PMV® 007 (Aqua color™)

27

---

---

---

---


---

---

---

### Benefits of a *Closed System*

- Restores airflow and sensation
  - Voice and communication
  - Sense of taste and smell
  - Eliminates need for finger occlusion or chin dropping



28

---

---

---

---

---

---

---

---

### Benefits of a *Closed System*

Impacts swallowing and may reduce aspiration

- Suck-Swallow-Breathe
- Airflow
- Sensation
- Positive airway pressure



29

---

---

---

---

---

---

---

---

### Benefits of a *Closed System*

- Restores physiologic positive pressure
  - Improved gas exchange
  - Improved oxygen saturation levels
  - Decreased risk of atelectasis



30

---

---

---

---

---

---

---

---

### Noah – Improved Cough



31

---

---

---

---

---

---

---

### Improved Secretion Management



1. Abraham, 2009

- Research study:
  - 24/49 children wearing PMV full time
  - Secretion management within normal limits in average of 2 weeks
  - Parent report QOL increases
  - Used valve for occasional sounds and secretion control

32

---

---

---

---


---

---

---

### Benefits of a *Closed System*

- Improves quality of life



33

---

---

---

---

---

---

---

### Pre-Assessment Considerations

- Inpatient vs Outpatient vs Home Health
- Factors to consider:
  - Parental involvement
  - Home environment
  - Proximity to therapy
  - Socioeconomic status
  - Insurance coverage



34

---

---

---

---

---

---

---

### Patient Selection Guidelines

- Awake and alert
- Medically stable
- Manage complete cuff deflation
- Manageable secretions
- Patent upper airway



35

---

---

---


---

---

---

---

### Application Checklist:



36

---

---

---


---

---

---

---

### Checklist: Educate



37

---

---

---

---

---

---

---

---

### Pediatrics

#### Facilitate Success: Play, Play, Play!

- Best in play-like environment (bubbles, singing, whistles, and more)
- Environment of trust
- Use distractions to your advantage
- Positive feedback



CHILD LIFE SPECIALIST, SIBLINGS,  
RECREATION THERAPISTS

38

---

---

---

---


---

---

---

---

### Warning Label



39

---

---

---

---

---

---

---

---

## Take Baseline Measurements

- Oxygenation
- Vital signs
- Breath sounds
- Color
- Work of breathing
- Patient responsiveness



40

---

---

---

---

---


---

---

---

## Patient Readiness: Adults and Peds

- Positioning
- State regulation
  - Clear sleep state
- Modulation
  - Robust crying
  - Periods of quiet alertness
- Attentional interaction
  - Attends to visual and auditory stimuli



41

---

---

---

---

---

---

---

---

## Suctioning

- Competencies
- When to suction
- How often
- Oral and tracheal
- Secretions?
  - Color
  - Smell
  - Thickness



42

---

---

---

---

---

---

---

---

## Assess Upper Airway Patency

- Deflate cuff
- Ask patient to inhale
- Finger occlude and speak or cough on exhalation



43

---

---

---

---


---

---

---

---

## Deflate Cuff



44

---

---

---

---

---

---

---

---

## Airway Patency for Valve Placement

### Adults

- Finger Occlusion
  - Describe breathing
    - Respiratory distress?
    - Work of breathing?
- TTP?

### Pediatrics

- Finger occlusion trials
  - Compliance
  - Tolerance
- TTP



45

---

---

---

---

---

---

---


---



## Assessment for Placement

### • Transtracheal pressures measurements

- Back pressure
- Air trapping
- Assessing for patent upper airway



46

---

---

---


---

---

---

---

## Candidate for Placement? Place PMV



47

---

---

---

---

---

---

---

## Assessment

### • What to watch for:

- Signs of Stability and Stress
  - Autonomic
    - Smooth breathing
    - Stable color
  - Motor
    - Good posture
    - Good tone



48

---

---

---

---


---

---

---

### After Placing the Valve

- Allow patient to adjust to airflow change
- Continue education and reassurance
- Establish phonation
- Continue to monitor for any changes from baseline measurements
- Remove valve if any significant changes occur



49

---

---

---

---

---

---

---

---

### Monitor Baseline Measurements

- Oxygenation
- Vital signs
- Breath sounds
- Color
- Work of breathing
- Patient responsiveness



50

---

---

---

---

---


---

---

---

### Following Placement

- Wear times vary - monitor
- Confer with medical team as needed
- Post warning labels



51

---

---

---

---

---


---

---

---

### Transitioning & Troubleshooting

- Excessive coughing
- Air trapping/back pressure
- Need for retraining
- Laryngeal/pharyngeal muscle rehab
- Psychological issues



52

---

---

---


---

---

---

---

### Air Whoosh – Back Pressure



53

---

---

---

---

---

---

---

### Factors Affecting Expiratory Air Flow

- Size or type of tracheostomy tube
- Presence and degree of obstruction
- Edema
- Secretions
- Incomplete cuff deflation
- Foam filled cuff
- Tube position



54

---

---

---

---

---

---

---

### Troubleshooting: Downsize or Different Brand



55

---

---

---


---

---

---

---

### Subglottic Stenosis



56

---

---

---

---

---

---

---

### Closed Tracheostomy Tube



57

---

---

---

---

---

---

---



58

---

---

---

---

---

---

---

#### Humidification With PassyMuir Valve

- Humidification is recommended
- Use with heat-moisture exchanger (HME) is ineffective
- Consider using with heated humidification
- Remove Valve for medicated aerosol treatments

59

---

---

---

---


---

---

---

#### Care and Cleaning

- Average lifetime of 2 months



60

---

---

---


---

---

---

---

### Supplemental Oxygen Adapter



PMA® 2000

61

---

---

---

---

---

---

---

### Impacts on Pressure and TTP Measurements



62

---

---

---

---

---

---

---

### Pressure: Resistance

Factor	Affected By
Length of the system	Constant; not a factor
Viscosity of air	Usually constant; humidity and altitude may alter viscosity slightly
Diameter of airways	Physical obstruction
Upper Airways	Mucous, stenosis, granulation tissue, and other factors
Lower Airways	Broncho constriction/dilation, tumors, disease, and other factors

63

---

---

---

---


---

---

---

### Primary Pressures Related to Airway Disorders

- Aerodigestive Tract
  - Respiratory
    - Thoracic cavity
    - Lungs
  - Digestive Tract
    - Oral
      - Intra-oral
      - Suck and swallow
    - Pharyngeal
      - Bolus transition
    - Esophagus
      - UES opening
      - Esophageal pressures
    - Stomach
    - Intestines
    - Abdominal cavity



64

---

---

---


---

---

---

---

### Toddler Anatomy



65

---

---

---


---

---

---

---

### Older Child Anatomy



66

---

---

---

---


---

---

---

## Pressurized System

- Restored or improved pressurized system:
  - Intraoral
  - Subglottic pressure
  - Respiratory – PEEP
  - Esophageal ??
  - Intrathoracic
  - Respiratory
- Leads to improved:
  - Feeding and swallowing
  - Cough and throat clear
  - Trunk support and postural control
  - Respiratory function



67

## Posture and Trunk Control

Primary pressure regulators


- Diaphragm
- Vocal folds/glottis (pressure regulator)
- Accessory muscles stabilize the spine and ribcage



Massery, M. (2006). Multisystem consequences of impaired breathing mechanics and/or postural control in cardiovascular and pulmonary physical therapy evidence and practice. (4th ed). In Frownfelter, D., & Dean, E., (Eds). St. Louis, MO: Elsevier Health Sciences, 695.

68

## Posture and Trunk Control



69



### Open Tracheostomy Tube



70

---

---

---


---

---

---

---

### Vocalization with Cuff Deflation: Leak Speech



71

---

---

---

---

---

---

---

### Above the Cuff Vocalization (ACV)



72

---

---

---

---

---

---

---

### Airway Patency: TTP Measurement



73

---

---

---

---


---

---

---

---

### Transtacheal Pressure Manometry: Decision Tree



74

---

---

---

---

---

---

---

---

### Transtacheal Pressure Measurement



75

---

---

---

---

---

---

---

---

### Pediatric Evaluation with TTP



76

---

---

---

---

---

---

---

### Breakout Sessions: Tracheostomy Tubes and PMVs, Cuff Management, and Mock Assessments



77

---

---

---

---

---

---

---

### Tracheostomy Tubes and PMVs



78

---

---

---

---

---

---

---

### Cuff Management



79

---

---

---

---

---

---

---

### Valve On and Off



80

---

---

---

---

---

---

---

### Mock Assessments



81

---

---

---

---

---

---

---

### Mock Assessment Patient Three

- 3 yr old male s/p MVC
- SCI – C 4-5
- Right lobe pneumothorax with chest tube
- Dysautonomia



82

---

---

---

---

---

---

---



**What additional information do you need?**

83

---

---

---

---

---

---

---



**You find out the pneumothorax is unresolved and has drainage. Could you use the PMV?**

84

---

---

---

---

---

---

---



**You find out the pneumothorax is primarily resolved and no drainage through the chest. So, you proceed with your evaluation. What would you do first?**

---

---

---

---

---

---

---

85



**How might dysautonomia impact your assessment?**

---

---

---

---

---

---

---

86

### Mock Assessment Patient One

- 16 yr old female with TBI s/p MVC
- Inpatient Rehab setting
- Self-extub x3
- Recovered pneumothorax
  
- Pt awake, alert, and mouthing words

---

---

---

---

---

---

---

87

Order for Assessment of  
Cognition and Speaking Valve



**Which would you conduct first?**

---

---

---

---

---

---

---

88



Patient demonstrates a patent airway and is  
able to use the PMV but has no voicing. What  
would you evaluate next?

---

---

---

---

---

---

---

89

### Mock Assessment Patient Two

- 44 yr old male s/p MVC with TBI
- Atelectasis
- Tracheoesophageal fistula

---

---

---

---

---

---

---

90



**What could be a cause of atelectasis? (mark all that apply)**

91

---

---

---

---

---

---

---

PMV and Communication



92

---

---

---

---

---

---

---



**What do you notice about the patient's voice/speech?**

93

---

---

---

---

---

---

---



LUNCH



& LEARN

---

---

---

---

---

---

---



# Lunch and Learn: Case Studies Questions and Discussion

---

---

---

---

---

---

---

1




---

---

---

---

---

---

---

2




---

---

---


---

---

---

---

3



4

---

---

---

---

---

---

---

---



5

---

---

---

---

---

---

---

---

Transtracheal Pressure  
and Behavioral Considerations



6

---

---

---

---

---

---

---

---

### Questions & Discussion

- Patent Upper Airway
- Behavioral Considerations
- Physiologic Considerations



7

---

---

---

---

---

---

---

Be back in 15 minutes



8

---

---

---

---

---

---

---

Treatment Considerations:  
Voicing, Swallowing, Mobility, and  
Impacts on Development



9

---

---

---

---

---

---

---

### Advantages of a Closed System vs Open Trach

- Open trach
  - Reduced airflow
  - Reduced positive airway pressure
- Closed Respiratory System
  - Allows graded exhalation and pressure regulation
    - Feeding and Swallowing
    - Posture and balance
    - Upper extremity force/strength



10

---

---

---

---

---

---

---

### Developmental Assessment Considerations

Preterm to  
4 months

- Assess oral motor functioning & feeding
- Vocalizations (evaluate for Valve use)
- Cooing and crying

4 to 12  
months

- Assess means to communicate
- Ability to vocalize
- Receptive language
- Parent-child interaction

11

---

---

---

---

---

---

---

### Developmental Assessment Considerations

8-24  
months

- Support communication systems and receptive language
- Motor coordination
- Ability to vocalize and cognition

24 months  
and older

- Progress in cognition
- Phonology, receptive/expressive language
- Change in communication system

Older  
Children

- Evaluate language, cognition, speech, and swallowing, as appropriate.
- Monitor respiratory function and valve use.

12

---

---

---

---

---

---

---

### Adult Considerations

#### Young Adult

- Support communication systems
- Evaluate speech-language and cognition as needed.
- Monitor respiratory function and Valve use.

#### Older Adult

- Support communication systems
- Evaluate speech-language and cognition as needed.
- Monitor respiratory function and Valve use.

13

---

---

---


---

---

---

---

### Development: Gross Motor & Oral Motor Skills



14

---

---

---

---

---

---

---

### Speech Effects in Pediatrics

- Subglottic and intraoral pressures
  - Bias Open vs Bias Closed
- Poor breath support
- Disruption of pre-linguistic language development
- Delayed onset of babbling stage
- Expressive and receptive language delays

Crucial period  
9 months/  
4.8 months

Hull et al. (2005); Roberts et al. (2019); Watters (2017); Zajac, D.J., Fornataro-Clerici, L., & Roop, T.A. (1999).

15

---

---


---

---

---

---

---



16

---

---

---

---

---

---

---

---

### Ability to Vocalize

- To feel and hear their own vocalizations
- Need opportunities to practice speech production skills

17

---

---

---

---

---

---

---

---

### Habilitate!

Early intervention:

Restores normal upper airway physiology as soon as possible

- Supports:
  - Speech and language development
  - Social interactions
  - Cognitive development

Jackson D and Albamonte S. (1994). Enhancing communication with the Passy-Muir valve. *Pediatric Nursing*. 20(2): 149-53.

18

---

---

---

---

---

---

---

---

## Facilitate Success and Co-Treatments

- Increase child's comfort with exhalation and encourage vocalization with activities such as:

- Microphone (echoes)
- Bubbles
- Horns or whistles
- Talking Tom App



- Valve wearing schedule

- Incorporate multidiscipline
- Optimize rehabilitation
  - Activities of daily living
  - Transfers
  - Exercise
  - Toileting

19

## Speech and Language Development

Avoid delays or decrease them with use of PMV


- Cooing
- Babbling



1. Simon, B.M., Fowler, S.M., and Handler, S.D. (1983) Communication development in young children with long-term tracheostomies: Preliminary report. *International Journal of Pediatric Otorhinolaryngology*, 6: 37-60.
2. Kadosh, K.W. and Stein, R.E. (1985) Chronic pediatric tracheostomy: Assessment and implications for habilitation of voice, speech and language in young children. *International Journal of Pediatric Otorhinolaryngology*, 9: 165-171.
3. Kertoy, M.K., Guest, C.M., Quast, E., and Luth-Lai, M. (1995). Speech and phonological characteristics of individual children with a history of tracheostomy. *Journal of Speech, Language and Hearing Research*, 42(3): 621-635.

20

## Speech and Language



21



What to do with a weak voice?



22

---

---

---

---


---

---

---

---

Increase Oral Exhalation:  
Retrain Upper Airway



- Techniques:
- Imitation
- Blowing
- Bubbles
- Whistles
- Horns, kazoos
- Pinwheels
- Straws
- Cotton balls

23

---

---

---

---


---

---

---

---

Swallow  
Physiology



24

---

---

---


---

---

---

---

---



25

---

---

---


---

---

---

---

---



26

---

---

---


---

---

---

---

---



27

---

---

---

---

---

---

---

---



28

---

---

---

---

---

---

---

---

### Postural Stability Study Findings:

"...balance strategies are likely to be disadvantaged if the ability to recruit glottal structures as part of dynamic postural control is compromised. This would be clinically meaningful for patients with a tracheostomy..."



Massery, M., Hagins, M., Stafford, R., Moerchen, V., Hodges, P.W. (2013). Effect of airway control by glottal structures on postural stability. *J Appl Physiol* (1985). 115(4):1483.

29

---

---

---

---

---

---

---

---

### Co-therapies: Mobility and Speech



30

---

---

---

---

---

---

---

---

### Mobility Considerations



31

---

---

---

---

---

---

---

### Respiratory Muscle Strength Training (RMST)

#### •Inspiratory Training Improves:

- Lung volumes
- Vocal fold opening

#### •Expiratory Training Improves:

- Cough strength
- Suprahyoid complex activation
- Vocal fold closure
- Breath support



32

---

---

---

---

---

---

---

### Functional Outcomes: What does the evidence show?

#### •Cough Effectiveness

- Improve inspiratory and expiratory pressures

#### •Swallow Effectiveness

- Expiratory muscle strength training is an effective intervention for impaired swallowing function in acute stroke patients with dysphagia



• Hutcherson, K.A., et al., (2017), Expiratory muscle strength training for radiation-associated aspiration after head and neck cancer: A case series. *Laryngoscope*, doi: 10.1002/lary.26845. [Epub ahead of print].

• Jang, H.A., Jin-Hwa, Ju, Young, S.W., Hee-Yang, C. and Kilhan, C. (2017). Effects of expiratory muscle strength training on swallowing function in acute stroke patients with dysphagia. *Journal of Physical Therapy Science*, 29(4), 609-612.

33

---

---

---

---


---

---

---

### Not Just for Adults - Functional Outcomes: Evidence?

- Risk of re-intubation
  - Neuromuscular weakness at time of extubation in children
- Improved functions
  - IMT improved: trunk control, respiratory muscle strength, ADLs, functional exercise capacity, and QOL in children with CP



• Khemani et al., 2017; Keles et al., 2018

34

---

---

---

---


---

---

---

---

### Therapy Goal Areas



- Valve use - waking hours with  $SpO_2 > 90\%$
- Swallowing & secretion management
- Re-establish intra-oral airflow management
- Speech and language development or rehabilitation
- Voicing
- Decannulation

35

---

---

---

---

---

---

---

---

### Toby Tracheapuppet™



36

---

---

---

---


---

---

---

---

## Pediatric and Adult Considerations for Ventilator Management



37

---

---

---

---

---

---

---

---

## Can Patients Drink While Ventilated?



38

---

---

---

---

---

---

---

---

## Peds vs Adults

### Adults = Trucks

- Stable, durable, and built for load
- Slower to change — takes time to see effects
- Can tolerate more variability
- More reserve — but slower recovery if overwhelmed
- Handle changes in fluid, meds, and settings with more cushion



39

---

---

---

---

---

---

---

---

## Peds vs Adults

### Pediatrics = Sports Cars

- High performance, low margin for error
- Quick to respond to changes — both good and bad
- Require precise, fine-tuned adjustments
- Less reserve — can tire out or crash quickly
- Small changes (fluid, meds, vent settings) can have big effects



40

---

---

---

---

---

---

---

## Indications for Invasive Mechanical Ventilation

- Can no longer support with NIV
- Airway protection
- Hypercapnic respiratory failure
- Hypoxemic respiratory failure
- Cardiovascular distress
- Anticipated patient decline or impending transfer



41

---

---

---

---

---


---

---

## Positive Pressure Ventilation

- The primary approach to mechanical ventilation is positive pressure
- The flow's driving pressure is higher than the pressure inside the lungs
- The pressure is high enough to overcome resistance

### Negative Pressure vs Positive Pressure Ventilation



42

---

---

---

---


---

---

---

## Invasive Ventilation

- Usually requires airway to be sealed with little to no leak present.
- Seal is achieved with a cuff at the end of the artificial airway.



43

---

---

---

---

---


---

---

---

## Cuff Deflation and Mechanical Ventilation

1. Set parameters do not change when cuff is deflated.
2. Cuff deflation generates less resistance to flow.
3. Ventilatory system is no longer sealed, there is a leak.



44

---

---

---

---

---

---

---

---

## Lung Mechanics

- Lung Compliance
- Factors affecting  $C_{LUNG}$ 
  - Pulmonary disease
  - Pulmonary edema or effusion
  - Tube position
  - Abdominal distension
- Airway Resistance
- Factors affecting  $R_{AW}$ 
  - Trach tube size
  - Secretions
  - Turbulent flow
  - High flow
  - Bronchospasm



45

---

---

---

---

---

---

---

---



### Why use a Passy-Muir Valve with patients who are mechanically ventilated?

- Verbal communication
- Improved lung recruitment and diaphragm involvement
- More rapid weaning from the ventilator
  - Rehabilitation tool
- Improved secretion management
  - More effective cough
  - Reduces need for suctioning
- Improves quality of life



46

---

---

---

---

---

---

---

---

### Can Patients Eat While Ventilated?



47

---

---

---

---

---

---

---

---

### Communication

- Adults:
  - Understand instructions & can express needs to report discomfort
- Pediatrics:
  - Non-Verbal Cues
    - Nasal Flaring
    - Grimacing
    - Smiling
  - Close monitoring



48

---

---

---


---

---

---

---

---



49

---

---

---


---

---

---

---

---



50

---

---

---


---

---

---

---

---



51

---

---

---

---

---


---

---

---

### Patient Selection Guidelines

- Awake and alert
- Medically stable
- Able to manage cuff deflation, if present
- Manageable secretions
- Low risk for gross aspiration
- Patent upper airway



52

---

---

---

---

---

---

---

---

### Conventional vs Non-Conventional Ventilation

#### Conventional Ventilation

- Follows a normal inhalation exhalation (normal breathing)
- Closely mirrors how one breathes without support

#### Non-Conventional Ventilation

- Does not follow normal breathing patterns
- Used when lungs are damaged – air is changed to protect them

53

---

---

---

---

---

---

---

---

### Conventional Modes of Ventilation

#### Control Modes

- Continuous Mandatory Ventilation (CMV)
- Assist Control (A/C)

#### Spontaneous Modes

- PS
- Continuous Positive Airway Pressure (CPAP)

#### Combination Mode

- Synchronized, Intermittent Mandatory Ventilation (SIMV)



54

---

---

---

---

---

---

---

---



### Modes of Ventilation

#### **Compatible with PMV Use:**

- A/C, VC, & PC
- SIMV
- CPAP/PS
- NIV

55

---

---

---

---

---

---

---



### Modes of Ventilation

#### **Not recommended with PMV:**

- PRVC
- APRV
- NAVA
- Other "auto" adjusted modes

56

---

---

---

---

---

---

---

### Ventilator Terms - Settings

- Mode of ventilation
- Tidal Volume
- Pressure
- Respiratory Rate
- Inspiratory Time
- Trigger sensitivity
- FiO<sub>2</sub>
- PEEP/CPAP



57

---

---


---

---

---

---

---



58

---

---

---

---

---


---

---

---

### Ventilator Terms – Alarm Settings

- High Pressure
- Low Pressure
- High and Low Exhaled Volumes
- High Respiratory Rate



59

---

---

---

---

---

---


---

---

### Assessing Ventilator Parameters

Three parameters that give you the general state of your patient's respiratory status:

1.  $\text{FiO}_2$
2. PEEP
3. PIP



60

---

---

---

---

---

---

---

---

## Step 2: Patient Preparation and Education

61

---

---

---

---

---

---

---

## Team Approach

- Timing and tube selection
- Introducing a speaking Valve
- When to downsize
- Plan of care
- Decannulation
- Impacts safety, length of stay, and cost



62

---

---

---

---

---

---

---



63

---

---

---

---

---

---

---

### Pre-Placement, General Observations, and Other Considerations

- Have a plan and block time
  - Pick a good time of the day
  - Reduce noise and interference
- Education
  - Reassure the patient
- Address pain issues
- Position the patient



64

---

---

---

---

---

---

---

---

### Pre-Placement, General Observations, and other Considerations

- Patient position
  - Congenital, anatomical differences, such as malformations, can affect the structure and function of the upper airway.
  - Proportionally larger heads
- Support with rolled towels or pillows



65

---

---

---

---

---

---

---

---

### Why Is Body Position Important?

- Tasks requiring trunk control:
  - Breathing
  - Coughing
  - Eating
  - Talking
  - Moving
  - Reaching
  - Bowel and bladder emptying



66

---

---

---

---

---

---

---

---

### Facilitate Success

- Environment of trust
  - Familiar voices
- Use distractions to your advantage
  - Soft toy
  - Pacifier
- Positive feedback
  - Time & trial



67

---

---

---

---

---

---

---

### Video



68

---

---

---

---

---

---

---

### Establish Baseline: Assess Vital Signs and Work of Breathing

- Oxygenation
- Vital Signs
- Breath sounds
- Color
- Work of breathing
- Patient responsiveness



69

---

---

---

---

---

---

---



### "Must Know" for PMV Use

- $\text{FiO}_2 \leq .50$
- $\text{PEEP} \leq 10 \text{ cmH}_2\text{O}$
- $\text{PIP} \leq 40 \text{ cmH}_2\text{O}$
- $\text{VTi}$  &  $\text{VTe}$
- Patient stability and ability to manage secretions and tolerate cuff deflation



70

---

---

---

---

---

---

---

---

### Assessing Ventilator Parameters



- $\text{FiO}_2$ 
  - Fraction of inspired Oxygen
  - Room Air 21%
  - Supplemental  $\text{O}_2 > 21\%$

71

---

---

---

---

---

---

---

---

### Assessing Ventilator Parameters



- PEEP
  - Positive End-Expiratory Pressure
  - Resistance to the exhaled volume that creates a back pressure that stents the alveoli open
  - PEEP and  $\text{FiO}_2$  work together to improve the oxygenation

72

---

---

---


---

---

---

---

---



### Assessing Ventilator Parameters

- PIP (Peak Inspiratory Pressure)
  - Max amount of pressure to deliver volume.
  - Sum of the inspiratory pressure required to deliver volume plus PEEP.
  - PIP indicates the compliance of the lungs.

73

---

---

---

---

---

---

---

---

### Step 3: Assess For Airway Patency

74

---

---

---

---

---


---

---

---

### Airway Patency Assessment with Mechanical Ventilation

- Requires complete cuff deflation
- Assess the leak or airflow into the upper airway
- Use vent parameters to determine airway patency
- The type of breath matters



75

---

---

---

---

---


---

---

---

## Cuff Deflation and Mechanical Ventilation

1. Set parameters do not change when cuff is deflated.
2. Cuff deflation generates less resistance to flow.
3. Ventilatory system is no longer sealed, there is a leak.



76

---

---

---

---

---

---

---

## Assessing Airflow Into the Upper Airway



77

---

---

---

---

---

---

---

## Upper Airway Assessment

- Transtracheal pressure (TTP) measurement



78

---

---

---

---


---

---

---

### Assessment Criteria

- Note PIP
- Observe inhaled and exhaled VT
- Achieve cuff deflation slowly, if present
- Assess volume changes (measure Vte)
- Listen for tracheal BS and assess leak through the upper airway



79

---

---

---

---

---

---

---

### Vent TTP Measurement



80

---

---

---

---

---

---

---

### TTP Measurement



81

---

---

---

---

---

---

---

## Q&A With Case Study Review and Troubleshooting: Pediatrics and Adults

82

---

---

---

---

---

---

---

### Case Study: Samuel

Home program  
focused on:

- Cooing
- Babbling
- Extended utterances
- Intelligibility
- Importance of positioning



83

---

---

---

---

---

---

---

### Case Studies: Low Level TBI

- Rehabilitation: Ranchos II-III
- Localized to generalized responses
- Stimulation



84

---

---

---

---

---

---

---

### Responsiveness



85

---

---

---

---

---

---

---

### Vocalization



86

---

---

---

---

---

---

---

### Case Study: Problem?



87

---

---

---

---

---

---

---

## Resolution



88

---

---

---

---

---

---

---

---

## Case Study: Grady



- Medical History
  - Born at 23 weeks
  - Periventricular Leukomalacia (PVL); Intraventricular Hemorrhage (IVH)
  - BPD
  - Bowel perforation
- DOL 155 – trach performed
- DOL 212 – PMV
- Vent settings: CPAP/PS
  - PEEP 8
  - PS 12
  - TTP >20

89

---

---

---

---

---

---

---

---

## Trial 1 – 15 minutes

- No vocalizations observed
- Patient calm, appeared happy and content
- Breathing unlabored with stable respiratory rate, no signs of distress
- Color pink, SpO<sub>2</sub> maintained within normal range
- Tolerance: good



90

---

---

---

---

---

---

---

---

### Trial 2 – 20 minutes

- Occasional cooing noted
- Patient alert, interacting with environment
- No signs of increased work of breathing
- Respiratory rate stable
- SpO<sub>2</sub> stable
- Tolerance: good



---

---

---

---

---

---

---

91

### Trial 3 – 15 minutes

- Increased vocalizations (cooing)
- Patient visually engaged and responsive
- Airflow through upper airway adequate
- RR unlabored
- V/S Stable
- Tolerance: good



---

---

---

---

---

---

---

92

### Trial 4 – 30 minutes

- Sustained vocalizations including crying and cooing
- First instance of audible crying
- Mother present, emotionally moved by hearing child's voice
- Patient demonstrated increased upper airway use
- No desaturation or signs of respiratory fatigue
- Tolerance: excellent



---

---

---

---

---

---

---

93



## Case Study: Grady



94

---

---

---

---

---

---

---

---

## What else should we consider?



- Manometry pressures (if applicable): for objective measurement of airway patency
- Secretion management: Was there increased secretion or need for suctioning?
- Upper airway patency: Was it visually assessed, or inferred from vocal quality?
- Cuff status (if trach is cuffed): e.g., cuff deflated
- Any changes in behavior or signs of fatigue post-trial

95

---

---

---

---

---

---

---

---

## Any questions?

- In summary, close the system to restore functions.



96

---

---

---

---

---

---

---

---



**Audience Q&A**

---

---

---

---

---

---

---

97

**Tomorrow:  
Ventilator Application of the Valve**



---

---

---

---


---

---

---

98

## Day Two Ventilator Application with the Valve



1

---

---

---

---

---


---

---

---

## Assessment Criteria

- Note PIP
- Observe inhaled and exhaled VT
- Achieve cuff deflation slowly, if present
- Assess volume changes (measure Vte)
- Listen for tracheal BS and assess leak through the upper airway



2

---

---

---

---

---

---

---

---

## "Must Know" for PMV Use

- $\text{FiO}_2 \leq .50$
- $\text{PEEP} \leq 10 \text{ cmH}_2\text{O}$
- $\text{PIP} \leq 40 \text{ cmH}_2\text{O}$
- VTi & VTe
- Patient stability and ability to manage secretions and tolerate cuff deflation



3

---

---

---


---

---


---

---

---



**RECOMMENDED**



**NOT RECOMMENDED**

**Modes of Ventilation**

**Compatible with PMV Use:**

- A/C, VC, & PC
- SIMV
- CPAP/PS
- NIV

4

---

---

---

---

---

---

---

---

VC: Patient Assessment


PIP 25 cmH<sub>2</sub>O

PEEP 5 cmH<sub>2</sub>O

FiO<sub>2</sub> 30%

VT 500 cc

Vte 500 cc



5

---

---

---

---

---


---

---

---

**Ventilator Assessment**

- Note Vent Settings:
  - Set Vt
  - PEEP
  - RR/f
  - FiO<sub>2</sub>
- Note Vent Measurements
  - PIP
  - Exhaled VT (Vte)
  - Total RR
  - MV



6

---

---

---

---

---

---

---

---

### Upper Airway Patency Assessment

- Turn Down PEEP
  - PEEP down by 5
- Then,
  - Slow cuff Deflation



7

---

---

---

---

---

---

---

### VC: Patient Assessment



PIP	25cmH <sub>2</sub> O
PEEP	8 cmH <sub>2</sub> O
FiO <sub>2</sub>	30%
VT	250cc
Vte	250cc

8

---

---

---

---

---

---

---

### PEEP

- Never turn the PEEP off
- Lower FiO<sub>2</sub> need
- Preserve FRC
- Prevent alveolar collapse
- Adjust by half



9

---

---

---

---


---

---

---

### Upper Airway Patency Assessment

- Turn Down PEEP
  - PEEP down by half-Peds
- Then,
  - Slow cuff Deflation



10

---

---

---

---

---

---


---

---

### Upper Airway Patency

Cuff Inflated-Closed Circuit

Cuff Deflated-Open Circuit



11

---

---

---

---


---

---

---

---

### PC: Patient Assessment



12

---

---

---

---

---

---

---

---

### Ventilator Assessment

- Note Vent Settings:

- Set IP
- PEEP
- RR/f
- $\text{FI}\text{O}_2$

- Note Vent Measurements

- PIP
- Exhaled VT (Vte)
- Total RR
- MV



13

---

---

---

---

---

---

---

---

### Upper Airway Patency Assessment

- Turn Down PEEP

- PEEP down by 5

- Then,
- Slow cuff Deflation



14

---

---

---

---

---

---

---

---

### PC: Patient Assessment



IP	20 cmH <sub>2</sub> O
PEEP	8 cmH <sub>2</sub> O
PIP	28
$\text{FI}\text{O}_2$	30%
Vte	250

15

---

---

---

---

---


---

---

---

### Upper Airway Patency Assessment

- Turn Down PEEP
  - PEEP down by half-Peds
- Then,
  - Slow cuff Deflation



16

---

---

---

---

---

---

---

---

### Ventilator Assessment

- Note Vent Settings:
  - Set IP
  - PEEP
  - RR/f
  - FIO<sub>2</sub>
- Note Vent Measurements
  - PIP
  - Exhaled VT (Vte)
  - Total RR
  - MV



17

---

---

---

---

---


---

---

---

### Upper Airway Patency Assessment

- Turn Down PEEP
  - PEEP down by half-Peds
- Then,
  - Slow cuff Deflation



18

---

---

---

---


---

---

---

---





19

---

---

---

---

---


---

---

---

### Distress May Appear As:

- Increased RR, HR
- Anxiety and fear
- Restlessness
- Increased irritability
- Stridor
- Decreased BS during auscultation
- Decreased chest movement



Squinting eyes

Flared Nostrils

Scrub Nail Beds

Tummy pulls in on back and chest

Rapid Belly Breathing

Pale Skin

20

---

---

---

---

---

---

---

---


### Troubleshooting

- PIP increases 2-4 cmH<sub>2</sub>O at a time
- Increased abdominal WOB noted

*What is wrong?*

- Remove the Valve and hear a "rush" of air coming from the trach

*What is this telling you?*



21

---

---

---

---

---

---

---

---

### Listen for Back Pressure



22

---

---

---

---

---

---

---



**What are some factors leading to back pressure?**

23

---

---

---

---


---

---

---

### Factors Affecting Upper Airway Patency

- Tube size and type
- Stenosis, edema, tumors, tracheomalacia
- Anatomical differences
- Foam cuff is an absolute contraindication.



24

---

---

---


---

---

---

---

### Tube Position is Important



25

---

---

---


---

---

---

---

### Tracheostomy Tubes and Impact on Airflow



26

---

---

---

---

---

---

---

### Step 4: Assemble the Necessary Parts & Pieces

27

---

---

---

---

---

---

---



**Which Valves can be used for in-line placement?**

---

---

---

---

---

---

---

28

### Ventilator Connections



---

---

---

---

---

---

---

29

### In-line Placement of the PMV® 2001 (Purple color™)



---

---

---


---

---

---

---

30



31

---

---

---

---

---

---

---

---



32

---

---

---


---

---

---

---

---



33

---

---

---

---

---

---

---

---

### Ventilator Adjustments

- Consider PEEP
- Evaluate Trigger Sensitivity
  - Pressure vs. Flow Trigger



34

---

---

---

---

---

---

---

### Ventilator Adjustments

- Volume Ventilation
- Increase VT in small increments to achieve pre-cuff deflation pressures (PIP)
  - Compensate for leak



35

---

---

---

---

---

---

---

### Pressure Adjustments

- Pressure ventilation
- If necessary, increase PC in small increments to achieve audible voice and adequate ventilation
  - Adjust I-Time



36

---

---

---

---


---

---

---

## Ventilator Assessment and Adjustments

- Pressure Support
  - Exp % sensitivity
  - Inspiratory cycle off
  - Set I-time
- Pressure Control
  - Set I-time



37

---

---

---

---

---

---

---

---

## Non-Invasive Ventilation(NIV)



- Approach to breath delivery is very similar to invasive ventilation.
- Difference is almost always have a leak present.

38

---

---

---

---

---

---

---

---

## Alarm Settings – Safe Practice

Low exhaled Vt and Ve alarms

Low pressure alarm


- Set 5 to 10 below PIP

High pressure alarm

- Set 5 to 10 above PIP

High respiratory rate

- 10 to 15 above baseline



39

---

---

---

---

---

---

---

---

### Humidification

- Heat/Moisture Exchanger (HME) is ineffective
- Use with Heated Systems
- Remove PMV for medicated treatments



40

---

---

---

---

---

---

---

### Case Study #1

- 8-month-old BPD
- Cuffless 3.5mm
- Awake and alert
- Hemodynamically stable



41

---

---

---

---

---

---

---

### Vent Settings

- SIMV 12
- IP 18 cmH<sub>2</sub>O
- PEEP 8 cmH<sub>2</sub>O
- PIP 26 cmH<sub>2</sub>O
- FIO<sub>2</sub> .24



42

---

---

---

---

---

---

---



### Airway Patency Assessment

- Transtracheal pressure measurement
- Observe several breaths
- Note pressure at the end of exhalation



43

---

---

---

---

---

---

---

---

### PassyMuir Valve Placement

- PEEP decreased to 4 cmH<sub>2</sub>O
  - Place PMV
  - PIP 22 cmH<sub>2</sub>O
  - Good cough and begins to make sounds
- Successful?



44

---

---

---

---

---

---

---

---

### Video



45

---

---

---

---

---

---

---

---

### Case Study #2

- Older child
- Trach 5.0, cuffed
- Awake and alert
- Hemodynamically stable



46

---

---

---

---

---

---

---

---

### Vent Settings

- SIMV 10
- VT 250 cc
- PIP 28 cmH<sub>2</sub>O
- PEEP 8 cmH<sub>2</sub>O
- FiO<sub>2</sub> .24
- PS 10 cmH<sub>2</sub>O



47

---

---

---

---

---

---

---

---

### Ventilator Assessment

- Note Vent Settings:
  - Set Vt
  - PEEP
  - RR/f
  - FiO<sub>2</sub>
- Note Vent Measurements
  - PIP
  - Exhaled Vt (Vte)
  - Total RR
  - MV



48

---

---

---

---

---


---

---

---

### Upper Airway Patency Assessment

- Turn Down PEEP
  - PEEP down by half-Peds
- Then,
  - Slow cuff Deflation



49

---

---

---

---

---

---

---

---

### Upper Airway Patency

Cuff Inflated-Closed Circuit



PEEP **8 cm H<sub>2</sub>O**  
 FiO<sub>2</sub> .24  
 PIP **28 cmH<sub>2</sub>O**  
 VT 250 cc  
 Vte **250 cc**

Cuff Deflated-Open Circuit



PEEP **4 cmH<sub>2</sub>O**  
 FiO<sub>2</sub> .24  
 PIP **23 cmH<sub>2</sub>O**  
 VT 250 cc  
 Vte **200 cc**

50

---

---

---

---

---

---

---

---

### Cuff Deflation Assessment

- Poor cough
- Not making any sounds

Would you place the Valve?



51

---

---

---

---

---

---

---

---

### Valve is NOT Placed In-Line

Why not?  
What are your  
recommendations?



52

---

---

---

---

---

---

---

### Re-assessment

Trach tube to size 4.5  
• Exhaled VT 100 cc  
• PIP 12 cmH<sub>2</sub>O

Would you place the valve?



53

---

---

---

---

---

---

---

### Ventilator Adjustments

- Patient is coughing
- Patient is verbalizing but doing so on inhalation and exhalation
- Appears to have difficulty cycling the breath
- Prolonged I-time

What would you adjust to help your patient?



54

---

---

---

---

---

---

---

### Troubleshooting

- Time or Flow Limit PS breath



55

---

---

---

---

---

---

---

### After Placement

- Treatment begins
  - Increase vocalizations
  - Secretion management
  - Address feeding and swallowing



56

---

---

---

---

---

---

---

### Review

- Assess ventilatory status
- Ensure airway patency
- Place Passy Muir Valve
- Ventilator Adjustment



57

---

---

---

---

---

---

---

## VITO Demonstration



58

---

---

---

---

---

---

---

---

## Hands-On: Parts and Pieces



59

---

---

---

---

---

---

---

---

## Hands-On: Inline Vent Placement



60

---

---

---

---

---

---

---

---

## Case Studies, Q&A, and Troubleshooting



61

---

---

---

---

---

---

---

## Sebastian – Case Study



62

---

---

---

---

---

---

---

## Charlie: 8 months old

- HRHS (dx in utero), white matter on brain, not eligible for transplant - d/c home hospice at 2.5 mos old
- Awake, alert, responsive, and active
- Cyanotic, SpO<sub>2</sub> 83–91% on high FIO<sub>2</sub>
- Tracheostomy present (3.0 uncuffed)
- LTV Vent settings:
  - PC = 12
  - PEEP = 4 cmH<sub>2</sub>O
  - SIMV = 25
- Spont. RR = 30/min, which is appropriate for age and consistent with the set SIMV rate.



63

---

---

---

---

---

---

---



### Concerns & Considerations Before PMV Placement

- **Cyanosis and variable SpO<sub>2</sub> (low 80s):**
  - poor oxygenation reserve.
- **Uncuffed tracheostomy tube (3.0):**
  - very small infants, tube size and airway resistance can become limiting.
- **Ventilator compatibility with PMV:**
  - The LTV ventilator can support PMV if appropriate settings are adjusted
- **HRHS and not transplant-eligible:**
  - Indicates limited long-term prognosis.
  - PMV may be used for quality of life, bonding, vocalization, or feeding coordination rather than long-term weaning.

64

---

---

---

---

---

---

---

---

### Pre-PMV Checklist

- Before attempting PMV, assess:
  - Upper airway patency-TTP
  - observe for distress, increased WOB, drop in SpO<sub>2</sub>
  - Team readiness (RT, speech-language pathologist, and provider) to monitor closely during initial trial.



65

---

---

---

---

---

---

---

---

### Off the vent

- **First Trial: Valve on the Vent**
  - TTP met the guidelines
  - Outpatient clinic setting
- **First Trial: Valve off the Vent**
  - Outpatient clinic and at home



66

---

---

---

---

---

---

---

---



## Case Study #2



67

---

---

---

---

---

---

---

**What would you check next?**

68

---

---

---

---

---

---

---

**What modes of ventilation  
are compatible with the  
Valve?**

69

---

---

---

---

---

---

---



**Cuff is deflated without adjusting PEEP. What are your considerations?**

70

---

---

---

---

---

---

---



**What is the first step in removing the Valve?**

71

---

---

---

---

---

---

---

Case Study #3



72

---

---


---

---

---

---

---



## Audience Q&A

73

---

---

---

---

---

---

---

## Thank you!

74

---

---


---

---


---

---

---



### Receiving CEUs for this Course



- You will have 5 days from the time this course ends to complete the evaluation, which is required to receive credit
- Go to: <https://ep.passy-muir.com>
- Login or create an account
- Click on the purple box
  - Upper righthand corner
  - Labeled "Enter Meeting Code Here"
- The meeting code is:

75

---

---

---

---

---

---

---