

Name:
MRN:
Date:

Speech and Language Therapy

Non-Ventilated Tracheostomy Weaning Plan



Cuff Deflation Regime



☐ Portex



☐ Shiley

Size _____

- ☐ Subglottic port
- ☐ XLT
- ☐ Cuffed
- ☐ Uncuffed
- ☐ Fenestrated
- ☐ Non- fenestrated

Time	Frequency
_____mins/hours	_____ times/day

Note _____

- Please monitor for ↓SpO₂, ↑RR, ↑Secretions → re-inflate cuff

One-Way Passy Muir Valve (PMV)



☐ PMV 007



☐ PMV 2001



☐ PMV 2000

Time	Frequency
_____mins/hours	_____ times/day

Note _____

- Please monitor for ↓SpO₂, ↑RR, ↑Secretions → remove PMV.
- **DO NOT** wear if cuff is inflated, on nebuliser or when sleeping.
- Clean the PMV daily according to manufacturer guidelines.

SLT:

Date:

Bleep: